

# Fluidigm® 48.770 Digital PCR Workflow Quick Reference

PN 68000147 Rev. B

For more information see, *BioMark Digital PCR Analysis Software User Guide, PN 68000100*


## 1 Priming the 48.770 Digital Array™ IFC



**CAUTION!** USE THE 48.770 CHIP WITHIN 24 HOURS OF OPENING THE PACKAGE.

- DUE TO DIFFERENT ACCUMULATOR VOLUMES, ONLY USE 48.48 SYRINGES WITH 300 µL OF CONTROL LINE FLUID.
- CONTROL LINE FLUID ON THE CHIP OR IN THE INLETS MAKES THE CHIP UNUSABLE.
- LOAD THE CHIP WITHIN 60 MINUTES OF PRIMING.

- 1 Inject control line fluid into each accumulator on the chip.
- 2 Remove and discard the blue protective film from the bottom of the chip.
- 3 Place the chip into the IFC (Integrated Fluidic Circuit) Controller MX, then run the **Prime (148x)** script to prime the control line fluid into the chip..



## 2 Preparing Sample Pre-Mix and Samples

- 1 Combine the components in the table below to make the Sample Pre-Mix and the final Sample Mixture (scale up appropriately for multiple runs).

SAMPLE PRE-MIX

Component	Volume per Inlet (µL)	Volume per Inlet with Overage (µL)	Volume per Chip (µL) (enough for 60 reactions)
TaqMan® Gene Expression Master Mix (Applied Biosystems, PN 4369016)*	2	3	180
20X GE Sample Loading Reagent (Fluidigm, PN 85000746) ●	0.4	0.6	36
20X gene-specific assays	0.2	0.3	18**
DNA-free water	0.2	0.3	18
DNA	1.2	1.8	
<b>Total</b>	<b>4</b>	<b>6</b>	<b>252</b>

\* TaqMan® Universal PCR Master Mix (Applied Biosystems, PN 4304437) may be substituted. Fluidigm recommends using TaqMan Gene Expression Master Mix for the Digital Array IFC.  
\*\* The 20X assay can be removed from the Sample Pre-Mix and added separately if different assays are to be used on the same chip.

- 2 In a DNA-free hood, combine the Taqman Gene Expression Master Mix, GE Sample Loading Reagent, DNA-free water and 20X assay(s) in a sterile tube—enough volume to fill the entire chip. 4.2 µL of this Sample Pre-Mix can then be aliquoted for each sample (48 total).
- 3 Remove these aliquots from the DNA-free hood and add 1.8 µL of DNA to each, making a total volume of 6 µL in each aliquot.

**NOTE** For a Copy Number Variation application, substitute RNase P for the DNA-free water.

## 3 Loading the Chip



**IMPORTANT!** VORTEX THOROUGHLY AND CENTRIFUGE ALL SAMPLE SOLUTIONS BEFORE PIPETTING INTO THE CHIP INLETS. FAILURE TO DO SO MAY RESULT IN A DECREASE IN DATA QUALITY.



**CAUTION!** WHILE PIPETTING, DO NOT GO PAST THE FIRST STOP ON THE PIPETTE. DOING SO MAY INTRODUCE AIR BUBBLES INTO THE INLETS.

- 1 When the **Prime (148x)** script has finished, remove the primed chip from the IFC Controller MX.
- 2 Pipette 10 µL of 1X GE Sample Loading Reagent into all hydration inlets.
- 3 Pipette 4 µL sample mix into the **sample inlets** on the chip.
- 4 Return the chip to the IFC Controller MX.
- 5 Using the IFC Controller MX software, run the **Load (148x)** script to load the samples into the chip.
- 6 When the **Load (148x)** script is complete, remove the chip from the IFC Controller MX.
- 7 Remove any dust particles or debris from the chip surface.


You are now ready for your chip run.



**CAUTION!** START THE CHIP RUN WITHIN 4 HOURS OF LOADING THE SAMPLES.

## 48.770 Pipetting Map


- Add 10 µL of 1X GE Sample Loading Reagent
- Pipette samples into 1-48



**A1**

## 4 Using the Data Collection Software

**IMPORTANT!** BE SURE TO SELECT ALL PROBE TYPES PRESENT IN YOUR EXPERIMENT. DATA ARE NOT COLLECTED ON UNSPECIFIED PROBES.


- 1 Double-click the Data Collection Software icon on the desktop to launch the software.
- 2 Click **Start a New Run**.
- 3 Check the status bar to verify that the lamp and the camera are ready. Make sure both are green before proceeding.  

- 4 Place the loaded chip into the reader.
- 5 Click **Load**.
- 6 Verify chip barcode and chip type.
  - a Choose project settings (if applicable).
  - b Click **Next**.
- 7 Chip Run file:
  - a Choose **New** or **Predefined**.
  - b Choose a file location for data storage.
  - c Click **Next**.
- 8 Application, Reference, Probes:
  - a Select Application Type—**Digital PCR**.
  - b Select Passive Reference (**ROX**).
  - c Select Assay—**Single probe, Two probes, or More than two probes**.
  - d Select probe types.
  - e Click **Next**.
- 9 Click **Browse** to find the appropriate thermal protocol file—**dPCR Standard v1.pcl**.
- 10 Confirm **Auto Exposure** is selected.
- 11 Click **Next**.
- 12 Verify the chip run information.
- 13 Click **Start Run**.

**NOTE** TO RUN THIS PROTOCOL AS AN END-POINT AND USE THE FLUIDIGM STAND-ALONE THERMAL CYCLER OR THE FLUIDIGM FC1™ CYCLER, REFER TO THE *FLUIDIGM STAND-ALONE THERMAL CYCLER USAGE QUICK REFERENCE* (PN 68000111) OR THE *FLUIDIGM FC1 CYCLER USAGE QUICK REFERENCE* (PN 100-1250), RESPECTIVELY.

## 5 Using the Digital PCR Analysis Software

**IMPORTANT!** BE SURE TO CLICK **ANALYZE** EACH TIME YOU CHANGE A PARAMETER IN THE SOFTWARE.

- 1 Double-click the Digital PCR Analysis software icon on the desktop to launch the software.
- 2 Click **Open a Chip Run**.
- 3 Double-click a **chiprun.bml** file to open it in the software.
- 4 Click **Sample and Detector Setup** in the Chip Explorer pane.
- 5 Click **New** or **Import**.
- 6 Highlight the wells and then annotate them.
- 7 Click **Editor** in the Sample and Detector Setup pane.
- 8 Choose **Sample Type** from the drop-down menu in the Editor.
- 9 Type a name for the sample.
- 10 Choose **Detector Type** from the drop-down menu in the Editor.
- 11 Type a name for the detector.
- 12 Click **Update** to see the changes reflected in the highlighted wells.
- 13 Click **Panel Summary** in the Chip Explorer pane.
- 14 Click **Analyze** in the Task pane.
- 15 Click **Panel Summary** or **Panel Details**.
- 16 Choose a view from the drop-down menu:
  - Results Table
  - Image View
  - Heat Map View



## Technical Support

### TELEPHONE

Within the United States: 1-866-358-4354

Outside the United States: 1-650-266-6100

### EMAIL

techsupport@fluidigm.com

© Copyright Fluidigm Corporation. All rights reserved. Fluidigm, the Fluidigm logo, Digital Array, FC1, and BioMark are trademarks or registered trademarks of Fluidigm Corporation in the U.S. and other countries. All other marks are the sole property of their respective owners. Refer to the *Digital PCR Analysis Software User Guide* (PN 68000100) for Fluidigm Patent, Limited License Agreement and Disclaimer. For research use only. Part Number 68000147 Rev. B Fluidigm recommends that you only purchase TaqMan® dual-labeled probes and/or other licensed PCR assay reagents from authorized sources. If you have any questions regarding whether you have a license to use particular reagents in PCR systems, you should contact the appropriate licensor and obtain clarification and their permission if necessary. For example, certain probes and their use may be covered by one or more patents held by Applied Biosystems and/or Roche Molecular Systems, which may be contacted at the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department, Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.