

## DELTAgene™ Assays


### KEY BENEFITS

- High quality—sensitivity and linearity similar to probe-based assays
- MIQE\* compliant—assay primer sequences provided
- Cost savings—low startup and running costs


Fluidigm DELTAgene Assays are high-quality, rapid turnaround qPCR gene expression assays. They permit you to take full advantage of your BioMark™ HD System with minimal experiment setup time and validated protocols that provide quality results.

Flexible groups of biologically-related genes are available for your specific requirements; there is no need to use fixed content. Simply provide us with your genes of interest (via RefSeq IDs) and the target species. All assays are ready-to-use via Fluidigm standard protocols. Assays are offered with or without wet-lab testing. We recommend using the second generation EvaGreen® dsDNA binding dye.

$C_q$  difference between 1, 10 and 100 cells easily distinguished



Single  $T_m$  Peak



Example Data: 1, 10 and 100 cell, custom EvaGreen®  
Assay Linearity Data (triplicates shown)

### PRODUCT DETAILS

- Primers sufficient for >100 chips; supplied in a 96-well plate
- Amplicons designed to cross an intron wherever possible (avoid amplifying genomic DNA)
- Minimum of 48 assays are designed to any RefSeq including human, mouse and rat
- Single-cell gene expression protocols available
- Turnaround time is three weeks for bioinformatically tested assays and six weeks for wet-lab tested assays
- Will design custom panels/pathways upon request

\* Minimum Information for Publication of Quantitative Real-Time PCR Experiments

## WORKFLOW

1 Assay Design	2 Prime	3 Transfer	4 Load	5 Run	6 Analyze
Fluidigm DELTAgene Assays are bioinformatically tested. Wet testing is optional.	Prime the Dynamic Array IFC to close the interface valves, preventing premature mixing of samples and assays.	Pipette samples, premixed with master mix, into separate sample inlets and the primer sets into separate inlets on the IFC.	Place the Dynamic Array IFC on the IFC Controller, and load the assay components into reaction chambers. Assay components are automatically combined on chip.	Place the Dynamic Array IFC on the BioMark HD System for thermal cycling and fluorescence detection.	Use real-time qPCR Analysis software to view and to interact with amplification curves, color-coded heat maps, and C <sub>q</sub> data for the run.

## ORDERING INFORMATION

Product	P/N
DELTAgene Assays (wet tested)	ASY-GE-WET
DELTAgene Assays (not wet tested)	ASY-GE

Biologically grouped assays are available for selection to build panels of your choice. Visit [www.fluidigm.com](http://www.fluidigm.com) for more details.

## FLUIDIGM SOLUTIONS FOR GENETIC ANALYSIS

## DELTAgene™ Assays

Custom assays designed for the BioMark™ HD System

## Dynamic Array™ IFCs

Consumable IFCs for high-throughput gene expression analysis and SNP genotyping

## Digital Array™ IFCs

Consumable IFCs for digital PCR

## IFC Controller

Integrated hardware/software for loading IFCs

## FC1™ Cyclor

Integrated hardware/software for thermal cycling of IFCs

## Software Suite

Analysis software for gene expression analysis, and automatic incorporation of C<sub>q</sub>

© 2011 Fluidigm Corporation. All rights reserved. Fluidigm, the Fluidigm logo, BioMark, DELTAgene, Digital Array, Dynamic Array and FC1 are trademarks or registered trademarks of Fluidigm Corporation in the U.S. and/or other countries. All other trademarks are the property of their respective owners. Fluidigm recommends that you only purchase licensed PCR assay reagents from authorized sources.

FOR RESEARCH USE ONLY.

100-3035 05/2011

## Corporate Headquarters

7000 Shoreline Court, Suite 100  
South San Francisco, CA 94080 USA  
Toll-free: 1.866.FLUIDLINE | Fax: 650.871.7152  
[www.fluidigm.com](http://www.fluidigm.com)

## Sales

North America | +1 650.266.6170 | [info-us@fluidigm.com](mailto:info-us@fluidigm.com)  
Europe/EMEA | +33 1 60 92 42 40 | [info-europe@fluidigm.com](mailto:info-europe@fluidigm.com)  
Japan | +81 3 3555 2351 | [info-japan@fluidigm.com](mailto:info-japan@fluidigm.com)  
Asia | +1 206.724.8386 | [info-asia@fluidigm.com](mailto:info-asia@fluidigm.com)

**Fluidigm** 